



LEED Product Data Sheet

Enterprise Adaptable Work Station

Recycled Content	Component Material	Total claim 20% Post Consumer Recycled Content. Of that 20% the ratio is broken down as follows:		
	18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge	Post Consumer 18%	Pre Consumer 2%	
Recycled Content	Earth Resin <i>(If Specified)</i>	Post Consumer 10%	Pre Consumer 0%	
Regional Materials	Regional Materials Call Kewaunee	<i>Must be within 500 miles to comply (Factory,Project, and Material supplier)</i>		
Adhesives and Sealants	Adhesives and Sealants Only applied on Site	Product Name Smooth On	VOC Content 0% g/L	VOC Content Allowed 250 g/L
Low emitting Materials	Paints and Coatings Not Applicable for V3, but applicable for V4	Call Kewaunee		

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com