



## LEED Product Data Sheet

### Fume Hoods

<b>Recycled Content</b>	<b>Component Material</b>	<b>Total claim 20% Recycled Content.</b>	
	18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge	<b>Of that 20% the ratio is broken down as follows:</b>	
	<b>STAINLESS STEEL</b>	<b>Post Consumer</b>	<b>Pre Consumer</b>
		18%	2%
		<b>Stainless Steel is 90% recycled content. Of that 90% the ratio is broken down as follows:</b>	
		80%	20%

*For Fume hoods use 60% of the total weight of the hood as your Steel amount.*

**CALCULATION: Recycled Content Value = (% Post Consumer Recycled Content x Material Cost) + 0.5 ( % Pre Consumer or Post Industrial Recycled Content x Material Cost)**

<b>Recycled Content</b>	Earth Resin ( <i>If Specified</i> )	<b>Post Consumer</b>	<b>Pre Consumer</b>
		10%	0%
<b>Regional Materials</b>	<b>Regional Materials</b>		
	Call Kewaunee	<i>Must be within 500 miles to comply (Factory, Project, and Material supplier)</i>	
<b>Adhesives and Sealants</b>	<b>Adhesives and Sealants</b>	<b>Product Name</b>	<b>VOC Content</b>
	Only applied on Site	Smooth On	0% g/L
			<b>VOC Content Allowed</b>
			250 g/L
<b>Low emitting Materials</b>	<b>Paints and Coatings</b>		
	Not Applicable for V3, but applicable for V4	<b>Call Kewaunee</b>	

Contact Rodney LaBelle for additional information. 704-871-3271 or [rodneylabelle@kewaunee.com](mailto:rodneylabelle@kewaunee.com)