



LEED Product Data Sheet

Research Collection

| | | | | |
|---|---|--|---------------------|----------------------------|
| Recycled Content | Component Material | Total claim 20% Recycled Content. | | |
| | 18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge | Of that 20% the ratio is broken down as follows: | | |
| | | Post Consumer | Pre Consumer | |
| | | 18% | 2% | |
| | STAINLESS STEEL | Stainless Steel is 90% recycled content. Of that 90% the ratio is broken down as follows: | | |
| | | 80% | 20% | |
| CALCULATION: Recycled Content Value = (% Post Consumer Recycled Content x Material Cost) + 0.5 (% Pre Consumer or Post Industrial Recycled Content x Material Cost) | | | | |
| Recycled Content | Earth Resin (<i>If Specified</i>) | Post Consumer | Pre Consumer | |
| | | 10% | 0% | |
| Regional Materials | Regional Materials | | | |
| | Call Kewaunee | <i>Must be within 500 miles to comply (Factory, Project, and Material supplier)</i> | | |
| Adhesives and Sealants | Adhesives and Sealants | Product Name | VOC Content | VOC Content Allowed |
| | Only applied on Site | Smooth On | 0% g/L | 250 g/L |
| Low emitting Materials | Paints and Coatings | | | |
| | Not Applicable for V3, but applicable for V4 | Call Kewaunee | | |

Contact Rodney LaBelle for additional information.
704-871-3271 or rodneylabelle@kewaunee.com