



LEED Product Data Sheet

Signature Series Style 05

	Component Material	Post Consumer recycled content	Pre Consumer recycled content
Recycled Content	Hardwood Plywood and Solids	0%	0%
	Particle Board	0%	100%
		Total claim 25% Pre Consumer Recycled Content	
Recycled Content	Earth Resin (<i>If Specified</i>)	Post Consumer	Pre Consumer
		10%	0%
Regional Materials	Regional Materials		
	Call Kewaunee	<i>Must be within 500 miles to comply (Factory, Project, and Material supplier)</i>	
FSC/COC	FSC certified	FSC Claim	Kewaunee FSC Certificate
	Upon Request	FSC Mix 75%	RA-COC-002756
Adhesives and	Adhesives and Sealants	Product Name	VOC Content
	Only applied on Site	Smooth On	0% g/L
			VOC Content Allowed
			250 g/L
Composite Wood	NAF		
	Upon Request		
Low emitting	Paints and Coatings		
	Not Applicable for V3, but applicable for V4	Call Kewaunee	

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com