

Enterprise Workstation Support Frame

FE30M8402XX-00

Single-sided, 2 inch Post Support Frame

Note:

Left Post will accept either electrical fixtures or service fittings but not both.

Note:

Right Post will accept either electrical fixtures or service fittings but not both.

Services on Left Post

Pre-piped to top of post
(see below for Piping Options)

Indicate service(s) required

indicate service required

indicate service required


indicate service required

Electric on Left Post

Pre-wired as selected below

Indicate outlet device type

- 120VAC Duplex
- 120VAC GFCI Duplex



Services on Right Post

Pre-piped to top of post
(see below for Piping Options)

Indicate service(s) required

indicate service required

indicate service required

indicate service required

Electric on Right Post

Pre-wired as selected below

Indicate outlet device type

- 120VAC Duplex
- 120VAC GFCI Duplex

Power Options: (indicate options required)
Electrical Device Color (pick one only)
 Black Ivory _____

Service Chase Receptacles (pick one only)
 Two (2) 120 VAC, 20 amp, Duplex
 Two (2) 120 VAC, 20 amp, GFCI, Duplex
 None

Service Chase/Post Electric Wiring (pick one only)
 3ft 4ft 5ft 6ft 8ft (above top of post)
 Single Circuit pigtail, twist-lock plug
 Double Circuit pigtail, twist-lock plug
 None

Electrical Pigtail Location (pick one only)
 Top of Left Upright
 Top of Right Upright

Piping Options: (indicate options required)
Piping Option (pick one only)
 Pipe to top of post

Piping/Hose End Option (pick one only)
 3/8" IPS Threaded Female (standard)
 3/8" IPS Threaded Male
 Keyed Male Quick-Connect W-8111-K_ _____

Specify Extension Hoses Separately

Support Frames are tested and certified to UL standard 61010-1 when pre-wired with cataloged fixtures.

FE30M8402 <input type="checkbox"/> <input type="checkbox"/> -00 <small>LENGTH</small>	<small>COLOR</small>	<small>ITEM NO.</small>	<small>PROJECT NAME:</small>
			<small>QUANTITY</small>