


# Enterprise Workstation Support Frame **FE32M8406XX-00**

Double-sided, Twin Tube Support Frame with 33½" high Vertical Chase



**Power Options:** (indicate options required)

**Electrical Device Color** (pick one only)

Black  Ivory  \_\_\_\_\_

**Service Chase Receptacles** (pick one only)

Two (2) 120 VAC, 20 amp, Duplex (each side)

Two (2) 120 VAC, 20 amp, GFCI, Duplex (each side)

None

\_\_\_\_\_

**Service Chase/Post Electric Wiring** (pick one only)

3ft  4ft  5ft  6ft  8ft (above top of post)

Single Circuit pigtail, twist-lock plug

Double Circuit pigtail, twist-lock plug

None

\_\_\_\_\_

**Electrical Pigtail Location** (pick one only)

Top of Left Upright

Top of Right Upright

**Piping Options:** (indicate options required)

**Piping Option** (pick one only)

Pipe to top of post

Pipe to top of post with extra:  
 3ft  4ft  5ft  6ft  8ft of hose

\_\_\_\_\_

**Piping/Hose End Option** (pick one only)

3/8" IPS Threaded Female (standard)

3/8" IPS Threaded Male

Keyed Male Quick-Connect W-8111-K\_

\_\_\_\_\_

**Specify Extension Hoses Separately**

**Support Frames are tested and certified to UL standard 61010-1 when pre-wired with cataloged fixtures.**

PROJECT NAME:	<b>FE32M8406</b> <input type="checkbox"/> <input type="checkbox"/> <b>LENGTH</b> <input type="checkbox"/> <input type="checkbox"/> <b>-00</b>	<b>COLOR</b>	<b>ITEM NO.</b>
			<b>QUANTITY</b>