


Evolution Workstation

FV20MXXXXXX-XX

9" Column, Double-sided, Starter Module



Accent Panel Option: (pick one only)
 0 Plastic Laminate
 1 Powder Coated Steel

Foot Option: (pick one only)
 0 1 1/2" Diameter Glide
 1 4" Swivel Locking Caster
 2 4" Heavy Duty Leveling Caster

Work Top Material: (pick one only)
 P Plastic Laminate T Phenolic Resin
 E ESD Laminate R Kemresin
 - No Work Top

Extended Length Option: (pick one only)
 0 Standard Length B Extended on Both Ends
 L Extended on Left End R Extended on Right End

Left Column Piping Option (pick one only)
Length above top of 84" Column
 3ft 4ft 5ft 6ft 8ft

Piping/Hose End Option (pick one only)
 3/8" IPS Threaded Female (standard)
 3/8" IPS Threaded Male
 Keyed Male Quick-Connect W-8111-K_ _____

Right Column Piping Option (pick one only)
Length above top of 84" Column
 3ft 4ft 5ft 6ft 8ft


Piping/Hose End Option (pick one only)
 3/8" IPS Threaded Female (standard)
 3/8" IPS Threaded Male
 Keyed Male Quick-Connect W-8111-K_ _____

FV20M HEIGHT <input type="checkbox"/> <input type="checkbox"/> DEPTH <input type="checkbox"/> <input type="checkbox"/> LENGTH <input type="checkbox"/> <input type="checkbox"/> ACCENT <input type="checkbox"/> <input type="checkbox"/> FOOT <input type="checkbox"/> <input type="checkbox"/> - <input type="checkbox"/> <input type="checkbox"/>	WORKSTATION COLOR	ACCENT COLOR	TOP COLOR
	PROJECT NAME:	ITEM NO.	ROOM NO.

Evolution Workstation

FV20MXXXXXX-XX

9" Column, Double-sided, Starter Module



Electrical Device Color (pick one only)

Black Ivory Grey Red _____

all devices same color unless noted above

Left Column Electric Wiring (pick one only)

- Single Circuit, Straight plug (NEMA 5-20)
- Single Circuit, Twist-lock plug (NEMA L5-20)
- Double Circuit, Twist-lock plug (NEMA 14-20)
- None
- _____

Pigtail Length: (pick one only)

3ft 4ft 5ft 6ft 8ft _____

Exit Post: (pick one only)

at Top at Bottom

Right Column Electric Wiring (pick one only)

- Single Circuit, Straight plug (NEMA 5-20)
- Single Circuit, Twist-lock plug (NEMA L5-20)
- Double Circuit, Twist-lock plug (NEMA 14-20)
- None
- _____

Pigtail Length: (pick one only)

3ft 4ft 5ft 6ft 8ft _____

Exit Post: (pick one only)

at Top at Bottom