


# Evolution Workstation

# FV30MXXXXXX-XX

12" Column, Single-sided, Starter Module



**Services on Left Column**  
Pre-piped to top of column  
(see below for Piping Options)

**Service Index**  
 AIR    GAS    VAC  
 AR    NIT  

**Hose Type**   **Color**  
 S/S    PVC   \_\_\_\_\_

**Service Index**  
 AIR    GAS    VAC  
 AR    NIT  

**Hose Type**   **Color**  
 S/S    PVC   \_\_\_\_\_

**Service Index**  
 AIR    GAS    VAC  
 AR    NIT  

**Hose Type**   **Color**  
 S/S    PVC   \_\_\_\_\_

**Services on Left Column**  
Pre-piped to top of column  
(see below for Piping Options)

**Service Index**  
 AIR    GAS    VAC  
 AR    NIT  

**Hose Type**   **Color**  
 S/S    PVC   \_\_\_\_\_

**Service Index**  
 AIR    GAS    VAC  
 AR    NIT  

**Hose Type**   **Color**  
 S/S    PVC   \_\_\_\_\_

**Service Index**  
 AIR    GAS    VAC  
 AR    NIT  

**Hose Type**   **Color**  
 S/S    PVC   \_\_\_\_\_

**Service Fittings:**  
 W-0605-00 Ball Valve  
 W-0606-00 Needle Point  
 \_\_\_\_\_

**Service Fittings:**  
 W-0605-00 Ball Valve  
 W-0606-00 Needle Point  
 \_\_\_\_\_

**Accent Panel Option:** (pick one only)  
 0 Plastic Laminate  
 1 Powder Coated Steel

**Foot Option:** (pick one only)  
 0 1 1/2" Diameter Glide  
 1 4" Swivel Locking Caster  
 2 4" Heavy Duty Leveling Caster

**Work Top Material:** (pick one only)  
 P Plastic Laminate    T Phenolic Resin  
 E ESD Laminate    R Kemresin  
 - No Work Top

**Extended Length Option:** (pick one only)  
 0 Standard Length    B Extended on Both Ends  
 L Extended on Left End    R Extended on Right End

**Left Column Piping Option** (pick one only)  
**Length above to of 84" Column**  
 3ft    4ft    5ft    6ft    8ft  
 \_\_\_\_\_

**Piping/Hose End Option** (pick one only)  
 3/8" IPS Threaded Female (standard)  
 3/8" IPS Threaded Male  
 Keyed Male Quick-Connect W-8111-K\_ \_\_\_\_\_  
 \_\_\_\_\_

**Right Column Piping Option** (pick one only)  
**Length above to of 84" Column**  
 3ft    4ft    5ft    6ft    8ft  
 \_\_\_\_\_


**Piping/Hose End Option** (pick one only)  
 3/8" IPS Threaded Female (standard)  
 3/8" IPS Threaded Male  
 Keyed Male Quick-Connect W-8111-K\_ \_\_\_\_\_  
 \_\_\_\_\_

<b>FV30M</b>	HEIGHT	DEPTH	LENGTH	ACCENT	FOOT	WORKSTATION COLOR	ACCENT COLOR	TOP COLOR
	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>
PROJECT NAME:						ITEM NO.	ROOM NO.	QUANTITY

# Evolution Workstation

# FV30MXXXXXX-XX

12" Column, Single-sided, Starter Module



**Electrical Device Color (pick one only)**  
 Black  Ivory  Grey  Red  \_\_\_\_\_  
*all devices same color unless noted above*

**Left Column Electric Wiring (pick one only)**  
 Single Circuit, Straight plug (NEMA 5-20)  
 Single Circuit, Twist-lock plug (NEMA L5-20)  
 Double Circuit, Twist-lock plug (NEMA 14-20)  
 None  
 \_\_\_\_\_

**Pigtail Length: (pick one only)**  
 3ft  4ft  5ft  6ft  8ft  \_\_\_\_\_

**Exit Post: (pick one only)**  
 at Top  at Bottom

**Right Column Electric Wiring (pick one only)**  
 Single Circuit, Straight plug (NEMA 5-20)  
 Single Circuit, Twist-lock plug (NEMA L5-20)  
 Double Circuit, Twist-lock plug (NEMA 14-20)  
 None  
 \_\_\_\_\_

**Pigtail Length: (pick one only)**  
 3ft  4ft  5ft  6ft  8ft  \_\_\_\_\_

**Exit Post: (pick one only)**  
 at Top  at Bottom