

# Supreme Air Venturi Fume Hood

# V05F28...

General Purpose Bench Hood with Vertical Rising Sash

(indicate fittings/fixtures required)  
(hole plugged if no fitting specified)

(indicate fittings/fixtures required)  
(hole plugged if no fitting specified)

[Option K - Fan/Blower Switch]

**Light/Sash Stop Controller**

[Option A1/A3 Alarm]

[Water]

[Water]

[Option P1 - Cord Port]

[Option P1 - Cord Port]

**120VAC GFCI Double Duplex**

**120VAC GFCI Double Duplex**

**Pre-piping:** (pick one if required)

- IN** Installed (no piping)
- PU** Pre-piped Up
- PD** Pre-piped Down
- SP** Pre-piped Single Point
  - Up - Left Side
  - Up - Right Side

**Receptacles:** (pick one)

- S** Spec. Grade GFCI
- H** Hospital Grade GFCI

**Receptacle Color:** (pick one)

- K** Black
- W** White
- R** Red
- V** Ivory
- G** Grey

**Wiring:** (pick if required)

- U** Pre-wired - Top of Hood

- Length:** (pick one only)  
 48"  60"  72"  96"
- Interior Depth:** (pick one only)  
 24"  30"  36"

**Required Options:** (pick one from each group)

- Interior Liner Material** (pick one)
- G** Kemglass (Fiberglass Reinforced Polyester)
  - K** KMER (Kewaunee Modified Epoxy Resin)
  - L** Type 316L Stainless Steel
  - S** Type 304L Stainless Steel
  - T** Phenolic Resin

- Sash Frame** (pick one)
- N** Frameless
  - M** Powder Coated Steel Frame
  - S** Type 304L Stainless Steel Frame

- Sash Glass** (pick one)
- G1** Laminated Safety Glass
  - G2** Tempered Safety Glass

- Service Fitting Valve Type** (pick one)
- F1** Rod Control Needle Valves
  - F3** Front Load Needle Valves

**Add-on Options:** (pick all required)

- A1** Air Alert 600 Alarm
- A3** Air Alert 300 Alarm (pick only one)
- C** Stainless Steel Duct Collar
- D** Distillation Rack - Preparation
- K** Fan/Blower Switch - (1hp motor rated)
- L** Sash Stop Label
- O** Stainless Steel Airfoil
- P1** Cord Port
- Q** Stainless Steel Sash Pulls
- R1** Auto Sash Return (pick only one)
- R3** Auto Sash Operator (pick only one)
- S** Safety Shield
- T** Tissue Screen

- V** Modified By-pass for VAV Control System

VAV Controller Manf: \_\_\_\_\_

Model: \_\_\_\_\_

Minimum CFM:

| Sash Height    | Inside Depth         | Overall Length       | Hood Liner           | Sash Frame           | Sash Glass           | Service Fittings     | Electrical Fixtures  | Add-on Options (separated by commas) | COLOR | ITEM NO. | QUANTITY |
|----------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|----------------------|--------------------------------------|-------|----------|----------|
| <b>V05F 28</b> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | <input type="text"/> | _____                                |       |          |          |