

Supreme Air Venturi Fume Hood

V06F28...

General Purpose Bench Hood with Combination Vertical Rising/Horizontal Sash

(indicate fittings/fixtures required)
(hole plugged if no fitting specified)

(indicate fittings/fixtures required)
(hole plugged if no fitting specified)

[Option K - Fan/Blower Switch]

Light/Sash Stop Controller

[Option A1/A3 Alarm]

[Water]

[Water]

[Option P1 - Cord Port]

[Option P1 - Cord Port]

120VAC GFCI Double Duplex

120VAC GFCI Double Duplex

Pre-piping: (pick one if required)

- IN** Installed (no piping)
- PU** Pre-piped Up
- PD** Pre-piped Down
- SP** Pre-piped Single Point
 - Up - Left Side
 - Up - Right Side

Receptacles: (pick one)

- S** Spec. Grade GFCI
- H** Hospital Grade GFCI

Receptacle Color: (pick one)

- K** Black
- W** White
- R** Red
- V** Ivory
- G** Grey

Wiring: (pick if required)

- U** Pre-wired - Top of Hood

Length: (pick one only)
 48" 60" 72" 96"
Interior Depth: (pick one only)
 24" 30" 36"

Required Options: (pick one from each group)

- Interior Liner Material** (pick one)
- G** Kemglass (Fiberglass Reinforced Polyester)
 - K** KMER (Kewaunee Modified Epoxy Resin)
 - L** Type 316L Stainless Steel
 - S** Type 304L Stainless Steel
 - T** Phenolic Resin

- Sash Frame** (pick one)
- M** Powder Coated Steel Frame
 - S** Type 304L Stainless Steel Frame

- Sash Glass** (pick one)
- G1** Laminated Safety Glass
 - G2** Tempered Safety Glass

- Service Fitting Valve Type** (pick one)
- F1** Rod Control Needle Valves
 - F3** Front Load Needle Valves

Add-on Options: (pick all required)

- A2** Air Alert 600 Alarm
- A3** Air Alert 300 Alarm (pick only one)
- C** Stainless Steel Duct Collar
- D** Distillation Rack - Preparation
- K** Fan/Blower Switch - (1hp motor rated)
- L** Sash Stop Label
- O** Stainless Steel Airfoil
- P1** Cord Port
- Q** Stainless Steel Sash Pulls
- R1** Auto Sash Return (pick only one)
- R3** Auto Sash Operator (pick only one)
- T** Tissue Screen

- V** Modified By-pass for VAV Control System

VAV Controller Manf: _____

Model: _____

Minimum CFM:

Sash Height	Inside Depth	Overall Length	Hood Liner	Sash Frame	Sash Glass	Service Fittings	Electrical Fixtures	Add-on Options (separated by commas)	COLOR	ITEM NO.	QUANTITY
V06F 28	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	_____			