


Supreme Air Venturi Fume Hood

ADA Bench Hood with Vertical Rising Sash

V15F28...

(indicate fittings/fixtures required)
(hole plugged if no fitting specified)

(indicate fittings/fixtures required)
(hole plugged if no fitting specified)



[Option K - Fan/Blower Switch]

[Option A1/A3 Alarm]

[Water]

[Water]

[Option P1 - Cord Port]

[Option P1 - Cord Port]

Light/Sash Stop Controller

120VAC GFCI Double Duplex

Pre-piping: (pick one if required)

Receptacles: (pick one)

- IN** Installed (no piping)
- PU** Pre-piped Up
- PD** Pre-piped Down
- SP** Pre-piped Single Point
 - Up - Left Side
 - Up - Right Side

- S** Spec. Grade GFCI
- H** Hospital Grade GFCI

Receptacle Color: (pick one)

- K** Black
- W** White
- R** Red
- V** Ivory
- G** Grey

Wiring: (pick if required)

- U** Pre-wired - Top of Hood

Length: (pick one only)
 48" 60" 72" 96"
Interior Depth: (pick one only)
 24" 30" 36"

Required Options: (pick one from each group)

Add-on Options: (pick all required)

Interior Liner Material (pick one)

- G** Kemglass (Fiberglass Reinforced Polyester)
- K** KMER (Kewaunee Modified Epoxy Resin)
- L** Type 316L Stainless Steel
- S** Type 304L Stainless Steel
- T** Phenolic Resin

Sash Frame (pick one)

- N** Frameless
- M** Powder Coated Steel Frame
- S** Type 304L Stainless Steel Frame

Sash Glass (pick one)

- G1** Laminated Safety Glass
- G2** Tempered Safety Glass

Service Fitting Valve Type (pick one)

- F2** Rod Control Ball Valves

- A1** Air Alert 600 Alarm
- A3** Air Alert 300 Alarm (pick only one)
- C** Stainless Steel Duct Collar
- D** Distillation Rack - Preparation
- K** Fan/Blower Switch - (1hp motor rated)
- L** Sash Stop Label
- O** Stainless Steel Airfoil
- P1** Cord Port
- Q** Stainless Steel Sash Pulls
- R1** Auto Sash Return (pick only one)
- R3** Auto Sash Operator (pick only one)
- S** Safety Shield
- T** Tissue Screen

V Modified By-pass for VAV Control System

VAV Controller Manf: _____

Model: _____

Minimum CFM:

Sash Height	Inside Depth	Overall Length	Hood Liner	Sash Frame	Sash Glass	Service Fittings	Electrical Fixtures	Add-on Options (separated by commas)	COLOR	ITEM NO.	QUANTITY
V15F28	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>				