


Supreme Air Venturi Fume Hood

V26F35...

LX Series 60" Interior Height Bench Hood with 35" Combination Vertical Rising/Horizontal Sash



Required Options: (pick one from each group)

Interior Liner Material (pick one)

- G** Kemglass (Fiberglass Reinforced Polyester)
- K** KMER (Kewaunee Modified Epoxy Resin)
- L** Type 316L Stainless Steel
- S** Type 304L Stainless Steel
- T** Phenolic Resin

Sash Frame (pick one)

- M** Powder Coated Steel Frame
- S** Type 304L Stainless Steel Frame

Sash Glass (pick one)

- G1** Laminated Safety Glass
- G2** Tempered Safety Glass

Service Fitting Valve Type (pick one)

- F1** Rod Control Needle Valves
- F3** Front Load Needle Valves

Add-on Options: (pick all required)

- A2** Air Alert 600 Alarm
- A3** Air Alert 300 Alarm (pick only one)
- C** Stainless Steel Duct Collar
- D** Distillation Rack - Preparation
- K** Fan/Blower Switch - (1hp motor rated)
- L** Sash Stop Label
- O** Stainless Steel Airfoil
- P1** Cord Port
- Q** Stainless Steel Sash Pulls
- R1** Auto Sash Return (pick only one)
- R3** Auto Sash Operator (pick only one)
- T** Tissue Screen

V Modified By-pass for VAV Control System

VAV Controller Manf: _____

Model: _____

Minimum CFM:

Sash Height	Inside Depth	Overall Length	Hood Liner	Sash Frame	Sash Glass	Service Fittings	Electrical Fixtures	Add-on Options (separated by commas)	COLOR	ITEM NO.	QUANTITY
V26F <input type="text" value="35"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	<input type="text"/>	_____			