




BasikBench[®]

SIMPLICITY at its finest



BasikBench®

- Offers the ideal starting point for a wide range of laboratory, cleanroom, and quality control activities
- Delivers a cost effective option for creating a functional workspace
- BasikBench is a value-priced 4-leg height-adjustable work bench



Configure your BasikBench®

Choose your bench, and whether you would like it to arrive knockdown or fully assembled. Then choose your accessories.

FEATURES

- Shelving adjustable in 1" increments
- Will allow for 36" mobile storage units
- Allows for up to four service fixtures
- Data ports standard on uprights
- SEFA 10, Category 3 – 1,000 pound capacity
- Top horizontal support channel for fixture hose to rest in
- Bolt together construction
- Ships knock down or assembled
- Worktop height adjustable on 1" increments

BasikBench® 36

Simplicity at its finest, BasikBench 4-leg laboratory bench comes in several sizes with choice of glides or leveling casters.



Shown above is an assembled table with top.

Table = BB363060-CA
Resin Top = BB36R013060-BK

C = Leveling Caster

G = Glide



Minimum Height

Maximum Height


27" with caster
25" with glide

40" with caster
38" with glide

Available in our standard metal paint colors.

12 Black	463 Iron Ore	242 Flint Grey	85 Satin Grey
94 True Blue	20 Sand Tan	61 Light Neutral	78 Snow White
103 Smokey Blue	90 Amazon	102 Agave	104 Sage
87 Mountain Blue	62 Driftwood	97 Invigorate	99 Rapture Red

Part Number	Depth	Length
BB362436-XX	24	36
BB362442-XX	24	42
BB362448-XX	24	48
BB362454-XX	24	54
BB362460-XX	24	60
BB362466-XX	24	66
BB362472-XX	24	72
BB363036-XX	30	36
BB363042-XX	30	42
BB363048-XX	30	48
BB363054-XX	30	54
BB363060-XX	30	60
BB363066-XX	30	66
BB363072-XX	30	72
BB363636-XX	36	36
BB363642-XX	36	42
BB363648-XX	36	48
BB363654-XX	36	54
BB363660-XX	36	60
BB363666-XX	36	66
BB363672-XX	36	72



C=CASTER
G=GLIDE
A= ASSEMBLED
K= KNOCKDOWN

BasikBench® 84



BasikBench® 84 expands above the worksurface with uprights. Allowing additional accesories and shelves to be mounted above the worksurface. Fixtures and data can also be added to each upright.

Shown above is an assembled table with top.

Table = BB843060-GA
Resin Top = BB84R013060-BK

C = Leveling Caster **G = Glide**



Minimum Height at Worksurface	Maximum Height at Worksurface
27" with caster	40" with caster
25" with glide	38" with glide
at top of Uprights	at top of Uprights
74.725" with caster	87.725" with caster
72.725" with glide	85.725" with glide

Available in our standard metal paint colors.

12 Black	463 Iron Ore	242 Flint Grey	85 Satin Grey
94 True Blue	20 Sand Tan	61 Light Neutral	78 Snow White
103 Smokey Blue	90 Amazon	102 Agave	104 Sage
87 Mountain Blue	62 Driftwood	97 Invigorate	99 Rapture Red

Part Number	Depth	Length
BB843036-XX	30	36
BB843042-XX	30	42
BB843048-XX	30	48
BB843054-XX	30	54
BB843060-XX	30	60
BB843066-XX	30	66
BB843072-XX	30	72
BB843636-XX	36	36
BB843642-XX	36	42
BB843648-XX	36	48
BB843654-XX	36	54
BB843660-XX	36	60
BB843666-XX	36	66
BB843672-XX	36	72


BB84 -

C=CASTER
G=GLIDE
A=ASSEMBLED
K=KNOCKDOWN

Tops Materials

Kemresin® Epoxy Resin

Kemresin® is perfectly suited for use in the lab setting. With ultra-smooth non-porous surface appearance, composed of 1 inch thick molded epoxy resin, these tops come in four different colors. Colors include black, grey, putty and slate. With resistance to many acids, alkalis and solvents, this work top is also resistant to both heat and physical use. These Kemresin® tops can be used wherever the highest physical and chemical resistance properties are required.



Grey = GR Slate = SL Putty = PT Black = BK

Part Number	Depth	Length	Part Number	Depth	Length
BB36R012436	24	36	BB84R013036	30	36
BB36R012442	24	42	BB84R013042	30	42
BB36R012448	24	48	BB84R013048	30	48
BB36R012454	24	54	BB84R013054	30	54
BB36R012460	24	60	BB84R013060	30	60
BB36R012466	24	66	BB84R013066	30	66
BB36R013072	30	72	BB84R013072	30	72
BB36R013036	30	36	BB84R013636	36	36
BB36R013042	30	42	BB84R013642	36	42
BB36R013048	30	48	BB84R013648	36	48
BB36R013054	30	54	BB8401R3654	36	54
BB36R013060	30	60	BB84R013660	36	60
BB36R013066	30	66	BB84R013666	36	66
BB36R013072	30	72	BB84R013672	36	72
BB36R013636	36	36			
BB36R013642	36	42			
BB36R013648	36	48			
BB36R013654	36	54			
BB36R013660	36	60			
BB36R013666	36	66			
BB36R013672	36	72			

BB36R -

BB84R -

Accessories

UPPER SHELF

Part Number	Depth	Length
BBUS1236	12	36
BBUS1242	12	42
BBUS1248	12	48
BBUS1254	12	54
BBUS1260	12	60
BBUS1266	12	66
BBUS1272	12	72
BBUS1536	15	36
BBUS1542	15	42
BBUS1548	15	48
BBUS1554	15	54
BBUS1560	15	60
BBUS1566	15	66
BBUS1572	15	72



This shelf can be installed to meet several heights via varied slots in the uprights of this bench.

LOWER SHELF

Part Number	Depth	Length
BBL0936	9	36
BBL0942	9	42
BBL0948	9	48
BBL0954	9	54
BBL0960	9	60
BBL0966	9	66
BBL0972	9	72
BBL1236	12	36
BBL1242	12	42
BBL1248	12	48
BBL1254	12	54
BBL1260	12	60
BBL1266	12	66
BBL1272	12	72
BBL1536	15	36
BBL1542	15	42
BBL1548	15	48
BBL1554	15	54
BBL1560	15	60
BBL1566	15	66
BBL1572	15	72



This lower shelf is installed across the bench's lower support.

CABLE MANAGEMENT


Part Number	Length
BBCM36	36
BBCM42	42
BBCM48	48
BBCM54	54
BBCM60	60
BBCM66	66
BBCM72	72



Cable management organizes power and data cables for an organized lab.

NOTE: ACCESSORIES ARE ORDERED AND SHIPPED SEPARATELY


MONITOR ARM



LCDARMBB

This monitor arm rotates 360 degrees at two points. With an adjustable head with universal monitor mount, this arm easily clamps onto any BasikBench.

CPU HOLDER



CPUHOLDER

The CPU holder attaches via a bracket mount and swivels side to side. The heavy duty holder is adjusted 12" to 20" high and 3.5" to 9" depth.

FIXTURES

Choose up to two fixtures per side. Fixtures can be installed on the assembled or knockdown version, as well as ordered separately. The AIR/VAC/NIT service fixtures include PVC hose at desired length. GAS service fixtures include stainless steel hoses at desired length. Tech sheet required for installed fixtures. Fixture and hose come preassembled terminated in a female swivel.

Hose Type	Part Number	Total Length
AIR/VAC/NIT	BBFXPV046-IN	46"
AIR/VAC/NIT	BBFXPV096-IN	96"
AIR/VAC/NIT	BBFXPV0168-IN	168"
GAS	BBFXSP046-IN	46"
GAS	BBFXSP096-IN	96"
GAS	BBFXSP0168-IN	168"

-IN only if installed

Length of hose available above the post can be found by subtracting 39 for top and 44 for bottom locations



Accessories

NOTE: ACCESSORIES ARE ORDERED AND SHIPPED SEPARATELY

HANGING RAIL FOR SUSPENDED CABINET


Part Number	Length
BBHR36	36
BBHR42	42
BBHR48	48
BBHR54	54
BBHR60	60
BBHR66	66
BBHR72	72



This hanging rail allows for the install of suspended cabinets.

BIN RAIL

Part Number	Length
BBBR36	36
BBBR42	42
BBBR48	48
BBBR54	54
BBBR60	60
BBBR66	66
BBBR72	72




This bin rail allows for easy installation of bins via a slotted rail. Bins can easily be moved and reorganized.

Note: Bin not included.


POWER BAR (20 AMP)

Available in Single or Dual circuit

Part Number	Length	Single or Dual Circuit (S/D)
BBPB2036S	36	S
BBPB2042S	42	S
BBPB2048S	48	S
BBPB2054S	54	S
BBPB2060S	60	S
BBPB2066S	66	S
BBPB2072S	72	S
BBPB2036D	36	D
BBPB2042D	42	D
BBPB2048D	48	D
BBPB2054D	54	D
BBPB2060D	60	D
BBPB2066D	66	D
BBPB2072D	72	D



Single Circuit (S)




Dual Circuit (D)

Power bars available in dual or single circuit and can be located on the uprights. Dual circuit has two cords and two power switches. The Single circuit has one cord and one power switch. Both come with twist lock plug.

ACCESSORIES

POWER STRIP MOUNTING RAIL

Part Number	Length
BBPSRL36	36
BBPSRL42	42
BBPSRL48	48
BBPSRL54	54
BBPSRL60	60
BBPSRL66	66
BBPSRL72	72




Power strip mounting rail allows for user supplied and installed power strip.

LED TASK LIGHTING

Part Number	Length
BBML20	20
BBML40	40

LED Task Lighting consists of a thin-profile, white, LED strip, 9.5 mm high by 33mm deep with a touch dimmer switch. It is UL listed and furnished with a 120 VAC to 24 VDC power supply with ten foot power cord. The light includes a magnetic mounting strip for mounting to steel shelving.

Part Number	Length	Fits under
X-030002-00	500mm - 19 ^{5/8}	29" - 40" Shelves
X-030003-00	1000mm - 39 ^{3/8}	46" - 47" Shelves



GLIDES AND CASTERS

Glides and casters come in sets of four.

Part Number
BBCAST
BBGLIDE

G = Glide



C = Leveling Caster



FLUSH WORK SURFACES



Storage Cabinets

Enough clearance for our 36" standing height, mobile, cabinets and accessory rail for hanging suspended cabinets. These cabinets are available in many different finishes as found in Kewaunee's Catalogs.



P.O. Box 1842 • Statesville, NC 28687-1842
Phone: (704) 873-7202 • Fax: (704) 873-5160
E-mail: kscmarketing@kewaunee.com
Web Site: www.kewaunee.com

Kewaunee Labway India Pvt. Ltd.
No. CA-9A, 2nd Floor, Jigani Link Road
Jigani, Bangalore 562 106, India
Phone: 9180-27826725 • Fax: 9180-27826724
Web Site: www.kewaunee.in

Kewaunee Labway Asia Pte. Ltd.
194 Pandan Loop, #06-22
Pantech Business Hub, Singapore 128 383
Phone: 65-6773-0288 • Fax: 65-6773-2322