



## LEED Product Data Sheet

### Alpha System Adaptable Furniture

**Recycled Content**      **Component Material**  
18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge

<i>Total Recycled Steel content ratio is broken down as follows: Steel weight makes up 60% of total Alpha weight</i>	
Post Consumer	Pre Consumer
30%	28.9%

**Recycled Content**      Aluminum

<i>Aluminum weight makes up 40% of total Alpha weight</i>	
Post Consumer	Pre Consumer
25%	50%

**Adhesives and Sealants**      **Regional Materials**  
Call Kewaunee

*Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)*

**Low emitting materials**      **Adhesives and Sealants**  
Only applied on Site

Product Name	VOC Content	VOC Content Allowed
Smooth On	0% g/L	250 g/L

**IEQ 4.2**      **Paints and Coatings**  
Not Applicable

Contact Rodney LaBelle for additional information. 704-871-3271 or [rodneylabelle@kewaunee.com](mailto:rodneylabelle@kewaunee.com)