



LEED Product Data Sheet

Enterprise Adaptable Work Station

Recycled Content
Component Material
 18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge

Total Recycled Steel content ratio is broken down as follows:

Post Consumer 30%	Pre Consumer 28.9%
--------------------------	---------------------------

Regional Materials
Regional Materials
 Call Kewaunee

Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)

Adhesives and Sealants
Adhesives and Sealants
 Only applied on Site

Product Name	VOC Content	VOC Content Allowed
Smooth On	0% g/L	250 g/L

Low emitting materials
Paints and Coatings
 Not Applicable

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com