



LEED Product Data Sheet

Fume Hoods

Recycled Content	Component Material	Total Recycled Steel content ratio is broken down as follows:	
	18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge	Post Consumer 30%	Pre Consumer 28.9%
	STAINLESS STEEL	Stainless Steel IF SPECIFIED is 82% recycled content. Of that 82% the ratio is broken down as follows: 63.77%18.29%	
For Fume hoods use 60% of the total weight of the hood as your Steel amount.			
CALCULATION: Recycled Content Value = (% Post Consumer Recycled Content x Material Cost) + 0.5 (% Pre Consumer or Post Industrial Recycled Content x Material Cost)			

Regional Materials	Regional Materials	<i>Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)</i>		
	Call Kewaunee			
Adhesives and Sealants	Adhesives and Sealants	Product Name	VOC Content	VOC Content Allowance
	Only applied on Site	Smooth On	0% g/L	250 g/L
Low emitting materials	Paints and Coatings			
	Not Applicable			

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com