



LEED Product Data Sheet

Research Collection

Recycled Content	Component Material	Total Recycled Steel content ratio is broken down as follows	
	16,18,20 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge	Post Consumer 30%	Pre Consumer 28.9%
	STAINLESS STEEL	Stainless Steel is 82% recycled content. Of that 82% the rate ratio is broken down as follows:	
		63.77%	18.29%

CALCULATION: Recycled Content Value = (% Post Consumer Recycled Content x Material Cost) + 0.5 (% Pre Consumer or Post Industrial Recycled Content x Material Cost)

Regional Materials **Regional Materials**
Call Kewaunee

Adhesives and Sealants	Adhesives and Sealants	Product Name	VOC Content	VOC Content Allowed
	Only applied on Site	Smooth On	0% g/L	250 g/L

Low emitting materials **Paints and Coatings**
Not Applicable

Contact Rodney LaBelle for additional information.
704-871-3271 or rodneylabelle@kewaunee.com

