



LEED Product Data Sheet

Research Collection

COMPONENT MATERIAL CABINET BODY
Recycled Content 16,18,20 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge

Total Recycled Steel content ratio is broken down as follows:

Post Consumer	Pre Consumer
30%	28.9%

STAINLESS STEEL CABINET BODY

Stainless Steel is 82% recycled content. Of that 82% the ratio is broken down as follows:

63.77%	18.29%
--------	--------

WOOD DOOR FRONTS

0%	0%
----	----

CALCULATION: $Recycled\ Content\ Value = (\% \text{ Post Consumer Recycled Content} \times Material\ Cost) + 0.5 (\% \text{ Pre Consumer or Post Industrial Recycled Content} \times Material\ Cost)$

Regional Materials **Regional Materials**
 Call Kewaunee

Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)

Adhesives and Sealants **Adhesives and Sealants**
 Only applied on Site

Product Name	VOC Content	VOC Content Allowed
Smooth On	0% g/L	250 g/L

Low emitting materials **Paints and Coatings**
 Not Applicable

Contact Rodney LaBelle for additional information.
 704-871-3271 or rodneylabelle@kewaunee.com