



LEED Product Data Sheet

Signature Series Style 05

	Component Material	Post Consumer recycled content	Pre Consumer recycled content
Recycled Content	Bamboo Plywood and Solids Particle Board	0% 0%	0% 100%
		Total claim 25% Pre Consumer Recycled Content	

Regional Materials	Regional Materials Call Kewaunee	<i>Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)</i>		
---------------------------	--	--	--	--

FSC/COC	FSC certified Upon Request	FSC Claim FSC Mix 75%	Kewaunee FSC Certificate NC-COC-002756	
----------------	--------------------------------------	---------------------------------	--	--

Adhesives and	Adhesives and Sealants Only applied on Site	Product Name Smooth On	VOC Content 0% g/L	VOC Content Allowed 250 g/L
----------------------	---	----------------------------------	------------------------------	---------------------------------------

Composite Wood	NAF Upon Request
-----------------------	----------------------------

Low emitting materials	Paints and Coatings Not Applicable
-------------------------------	--

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com