



LEED Product Data Sheet

Signature Series Style 05

	Component Material	Post Consumer recycled content	Pre Consumer recycled content
Recycled Content	Hardwood Plywood and Solids	0%	0%
	Combi Core Plywood Particle Board is typically 85% of total Cabinet weight	0%	100%

Regional Materials

Regional Materials

Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)

FSC/COC

Call Kewaunee

FSC Claim

FSC Mix 75%

Kewaunee FSC Certificate

NC-COC-002756

FSC certified

Upon Request

Adhesives and

Adhesives and Sealants Only

applied on Site

Product Name

Smooth On

VOC Content

0% g/L

VOC Content Allowed

250 g/L

Composite Wood

NAF

Upon Request

Low emitting materials

Paints and Coatings

Not Applicable

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com