



LEED Product Data Sheet

Venturi Fume Hoods

Recycled Content	Component Material	Total Recycled Steel content ratio is broken down as follows:	
	18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge	Post Consumer 30%	Pre Consumer 28.9%
	STAINLESS STEEL (IF SPECIFIED)	Stainless Steel is 82% recycled content. Of that 82% the ratio is broken down as follows:	
		63.77%	18.29%

For Venturi Fume hoods use 60% of the total weight of the hood as your Steel amount.

CALCULATION: Recycled Content Value = (% Post Consumer Recycled Content x Material Cost) + 0.5 (% Pre Consumer or Post Industrial Recycled Content x Material Cost)

Regional Materials	Regional Materials Call Kewaunee	<i>Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)</i>		
Adhesives and Sealants	Adhesives and Sealants Only applied on Site	Product Name Smooth On	VOC Content 0% g/L	VOC Content Allowed 250 g/L
Low emitting materials	Paints and Coatings Not Applicable			

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com