

Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame

FEXXC8402XX-0

Power Cord (Max 2 per post)

Power Cord (Max 2 per post)

Raceway Electric (Max 10 per bench)

Raceway Electric (Max 10 per bench)

Services Left Post

Pre-piped to Top of Post  
See Catalog for Part Numbers

Indicates Fittings/Fixtures Required  
*Plugged if not selected*

indicate service required

indicate service required

indicate service required

Services Left Post

Pre-piped to Top of Post  
See Catalog for Part Numbers

Indicates Fittings/Fixtures Required  
*Plugged if not selected*

indicate service required

indicate service required

indicate service required

Electric on Left Post

Electric on Left Post

FE \_\_ C8402 \_\_ -

PROJECT NAME:

COLOR


ITEM NO

QUANTITY

# Enterprise Workstation Support Frame

Blank, 2 inch Round or Square Post Support Frame

# FEXXC8402XX-0



## SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

-Provide hole patterns in each post for edison iii attachment  
-Replace standard levelers with stainless levelers  
-Add casters (requires counter-balanced fe80 to maintain ul)  
-Reposition the standard horizontal raceway  
-Add (1) additional horizontal raceway  
-Change the single horiz. Raceway to a double-channel horiz. Raceway  
-Frame height mod (add or reduce height below lowest service fixture)

-Extend frame length (longer than 72" w/96" max) - centerpost remains  
-Reduce frame length (less than 54" w/24" min)- centerpost removed  
-Receive a cutout when no cataloged electric or fixture is ordered  
-Add (1) leg electric in addition to those offered standard  
-Replace leg receptacle with a data device  
-Add a non-cataloged electric or data device - **VOIDS UL**  
-Stainless construction) - **VOIDS UL**

## TO BE DEFINED ON A "-0" TECH SHEET:

-Data, elec, fixture c/o positions  
-Data, elec. Cord exit  
-Number of circuits and positions

-Modified frame dimensions  
-Added horiz. Raceway  
-Repositioned horiz. Raceway

FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

0324