

Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame

FEXXC8406XX-0

Power Cord (Max 2 per post)

Raceway Electric (Max 10 per bench)

Services Left Post

Pre-piped to Top of Post  
See Catalog for Part Numbers

Indicates Fittings/Fixtures Required  
*Plugged if not selected*

*indicate service required*

*indicate service required*

*indicate service required*

Electric on Left Post

Power Cord (Max 2 per post)

Raceway Electric (Max 10 per bench)

Services Left Post

Pre-piped to Top of Post  
See Catalog for Part Numbers

Indicates Fittings/Fixtures Required  
*Plugged if not selected*

*indicate service required*

*indicate service required*

*indicate service required*


Electric on Left Post

FE __ C8406 __ -	COLOR	ITEM NO
PROJECT NAME:		QUANTITY

# Enterprise Workstation Support Frame

Blank, 6 inch Round or Square Post Support Frame

# FEXXC8406XX-0



SOME EXAMPLES OF MODIFICATIONS REQUIRING A ZZ: (Consult Special Pricing for other modifications)

- Provide hole patterns in each post for edison iii attachment
- Replace standard levelers with stainless levelers
- Add casters (requires counter-balanced fe80 to maintain ul)
- Reposition the standard horizontal raceway
- Add (1) additional horizontal raceway
- Change the single horiz. Raceway to a double-channel horiz. Raceway
- Modify frame to receive vacuubrand fixture

- Frame height modification
- Extend frame length (longer than 72" w/96" max) - centerpost remains
- Reduce frame length (less than 54" w/24" min)- centerpost removed
- Receive a cutout when no cataloged electric or fixture is ordered
- Add (1) leg electric in addition to those offered standard
- Replace leg receptacle with a data device
- Add a non-cataloged electric or data device - **VOIDS UL**
- Stainless construction) - **VOIDS UL**

TO BE DEFINED ON A "-0" TECH SHEET:

- |                                     |                              |
|-------------------------------------|------------------------------|
| -Data, elec., fixture c/o positions | -Modified frame dimensions   |
| -Data, elec. Cord exit              | -Added horiz. Raceway        |
| -Number of circuits and positions   | -Repositioned horiz. Raceway |

FOR ANY SPECIAL DESIGN REQUIREMENT THAT IS NOT DESCRIBED ON THIS PAGE, CONTACT YOUR DEALER SUPPORT MANAGER WHO WILL ARRANGE A DESIGN CONSULTATION DISCUSSION

0324