



LEED Product Data Sheet

Evolution System Adaptable Furniture

Recycled Content **Component Material**
18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge

Total Recycled Steel content ratio is broken down as follows:Steel weight makes up 60% of total Evolution weight	
Post Consumer 30%	Pre Consumer 28.9%

Recycled Content Aluminum

Aluminum weight makes up 40% of total Evolution weight	
Post Consumer 25%	Pre Consumer 50%

Adhesives and Sealants **Regional Materials**
Call Kewaunee

Must be within 500 miles V3 or 100 miles V4 to comply (Factory,Project, and Material supplier)

Low emitting materials **Adhesives and Sealants**
Only applied on Site

Product Name	VOC Content	VOC Content Allowed
Smooth On	0% g/L	250 g/L

IEQ 4.2 **Paints and Coatings**
Not Applicable

Contact Rodney LaBelle for additional information. 704-871-3271 or rodneylabelle@kewaunee.com