



## LEED Product Data Sheet

### Basik Bench Work Station

<b>Recycled Content</b>	<b>Component Material</b>	<b>Total Recycled Steel content ratio is broken down as follows:</b>		
	18 Gauge furniture grade sheet metal with structural reinforcement members of 12 and 14 gauge	<b>Post Consumer</b>		<b>Pre Consumer</b>
		2.18%		10.19%
<b>Regional Materials</b>	<b>Regional Materials</b>	<i>Must be within 500 miles V3 or 100 miles V4 to comply (Factory, Project, and Material supplier)</i>		
	Call Kewaunee			
<b>Adhesives and Sealants</b>	<b>Adhesives and Sealants</b>	<b>Product Name</b>	<b>VOC Content</b>	<b>VOC Content Allowed</b>
	Only applied on Site	Smooth On	0% g/L	250 g/L
<b>Low emitting materials</b>	<b>Paints and Coatings</b>			
	Not Applicable			

Contact Rodney LaBelle for additional information. 704-871-3271 or [rodneylabelle@kewaunee.com](mailto:rodneylabelle@kewaunee.com)